

BEST AVAILABLE COPY

mix 200 μ l aliquots to gether.

Work cells are w/ PSE

Add final media (complete w/ serum)
Add final media of the DNA c/a

I devised strategies for making the Adeno ~~EE~~ sp1A / IB constructs.
Basically - move in the reported gene and then PCR the
minimal enhancer and move in it a second cloning.

For example, PSE driving p-gal in Adeno

orientation w/ Clal

1. Take XbaI fragment (3.88 kb) from pCMV-pot-gal
which has SV40 3' / 5', p-gal + polyA into
XbaI vector pAE1 sp1A / IB

orientation w/ Clal

2. PCR the minimal enhancer w/ HindIII ends. Note
HindIII does not cut in the 3.88 kb of p-gal.
Clone in the enhancer & we are done

PSE driving CAT in Adeno

orientation w/ EcoRI

1. Take XbaI/BamHI of pCAT basic and move
into XbaI/BamHI of pAE1 sp1A / IB, (1.7 kb)

orientation w/ Clal \pm EcoRI

2. In the resultant plasmids clone in HindIII ends
off the minimal PSE. The same as above...
HindIII does not cut CAT & only cuts IX in
the MCS of the vectors

PSE driving DTA in Adeno

orientation w/ EcoRI

1. Take XbaI/BamHI of 3'44 (1.4 kb) and clone into
XbaI/BamHI of pAE1 sp1A / IB

orientation w/ Clal \pm EcoRI

2. Clone in minimal PSE w/ HindIII ends
as above.

To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date

Recorded by

D. Henderson