

# BEST AVAILABLE COPY

mix 200  $\mu$ l aliquots to gether.  
 Work cells in w/ PSE  
 Add final media (complete w/ serum)  
 Add final media of the DNA c/a

I devised strategies for making the Adeno ~~EE~~ sp1A / IB constructs.  
 Basically - move in the reported gene and then PCR the  
 minimal enhancer and move in as a second cloning.

For example, PSE driving p-gal in Adeno

orientation w/ Clal

1. Take XbaI fragment (3.88 kb) from pCMV-pot-gal which has SV40 3' / 5', p-gal + polyA into XbaI vector pAE1 sp1A / IB
2. PCR the minimal enhancer w/ HindIII ends. Note HindIII does not cut in the 3.88 kb of p-gal. Clone in the enhancer & we are done

orientation w/ Clal

PSE driving CAT in Adeno

orientation w/ EcoRI

1. Take XbaI/BamHI of pCAT basic and move into XbaI/BamHI of pAE1 sp1A / IB, (1.7 kb)
2. In the resultant plasmids clone in HindIII ends off the minimal PSE. The same as above... HindIII does not cut CAT & only cuts 1X in the MCS of the vectors

orientation w/ Clal  $\pm$  EcoRI

PSE driving DTA in Adeno

orientation w/ EcoRI

1. Take XbaI/BamHI of 3'44 (1.4 kb) and clone into XbaI/BamHI of pAE1 sp1A / IB
2. Clone in minimal PSE w/ HindIII ends as above.

orientation w/ Clal  $\pm$  EcoRI

To Page No. \_\_\_\_\_

Witnessed & Understood by me,

Date

Invented by

Date

Recorded by

D. Henderson